

Plug*

USBPW



10 port plug*

USB10X



20 port plug*

USB20X



Charging and download cable

CheckmeCB



***Visi-Download for Viatom**



Replacement strap for MAX

540131



Replacement metal clip for strap

540129



Checkme O2 MAX oximeter
CheckmeO2MAX





Ring SpO₂ sensor with vibration

Reusable, cable length 130mm, 540124



*** Adult soft SpO₂ sensor,**

cable length 145mm, 540123



*** Paediatric soft SpO₂ sensor**

Reusable, cable length 130mm, 540126



*** Single-patient foam SpO₂ sensor**

Single patient use, 29mm between sensor and detector, cable length 130mm, 540127



*** Neonatal non-woven SpO₂ sensor**

Single patient use, 20mm between sensor and detector, cable length 130mm, 540128

* OPTIONAL ITEM / ACCESSORY ** ALTERNATIVES AVAILABLE

